

Bionic Buttonhole Needle



The gentle puncture

Besides the known advantages of the Bionic Fistula Needles, for example the unique design of the Bionic hub, the silicone coating, the low flow resistance and the slit-shaped back eye the Bionic Buttonhole Needle offers the perfect balance between a pointed and a blunt needle. The Bionic Buttonhole Needle which is especially developed for the Buttonhole Technique is a reasonable alternative particularly for home dialysis patients and patients who do self-puncturing.


Perfect Balance – the Bionic Buttonhole Needle



The perfect balance between pointed and blunt needles developed especially for the Buttonhole Technique.



Rotating wings and coloured markers for the perfect puncturing.



Red and blue clamps are easy to handle and offer safety.

The advantages at a glance:

- Functional design of the needle hub
- Rotating wings and coloured markers for the exact bevel positioning
- Colour coded light-weight clamps for safe handling
- Luer-Lock connector with a rough inner surface
- Reduced number of failed punctures
- Less pain experience
- Shortened bleeding time

Bionic Buttonhole Needle Product Range

Dimensions Ø	Ø x Length	Wing Colour	Order Code
16 G	1.6 x 20 mm		646-AGAM-BH
16 G	1.6 x 20 mm		646-VGAM-BH
16 G	1.6 x 25 mm		656-AGAM-BH
16 G	1.6 x 25 mm		656-VGAM-BH
15 G	1.8 x 20 mm		846-AGAM-BH
15 G	1.8 x 20 mm		846-VGAM-BH
15 G	1.8 x 25 mm		856-AGAM-BH
15 G	1.8 x 25 mm		856-VGAM-BH
14 G	2.0 x 20 mm		046-AGAM-BH
14 G	2.0 x 20 mm		046-VGAM-BH
14 G	2.0 x 25 mm		056-AGAM-BH
14 G	2.0 x 25 mm		056-VGAM-BH

A = Arterial with back eye, red clamp V = Venous without back eye, blue clamp Tube length 30 cm

