

Bionic-Accessories


secure, safe and effective

- for a secure dialysis


Universal cap, circulation connector, infusion systems, Single-Needle-Adapter, scissor clamps

- for an effective connection and disconnection

Sets for single and double lumen catheters, customized sets on request




Bionic-Accessories



AV-EXT2

Extension lines for dual lumen catheters, each 10 cm long. The clamps in red and blue assure safe and fast clamping.



AV-EXT2-Y30

Connection aid for fistula needles, each 30 cm long. Arterial line with y-fitting to add saline solution at the same time.



SNA3

Y-Adapter with two connection tubes, each 10 cm long. The colour coded clamps assure fast and safe clamping.



CLAMP-SP-02

Stop&Protect - simple disposable compression strap, practical dispenser box with 50 pieces



CATH-TAS

DIACATH catheter dressing with pocket made of non-woven fleece for easy and quick fixation of catheters.



K2565

Plastic scissor clamp for easy clamping of tubes of up to 7 mm diameter.



KPA06006-4 / KPA06004-4

Connection and disconnection set for single or double lumen catheter, contains different components, e.g. drape.



KPA06021-6

Connection and disconnection set for shunt puncture; contains various components, e.g. cover cloth, plaster.



KAB75

Sharps container with a volume of 7.5 liters.



HEDclip

HEDclip for strain relief of all common blood tube systems with a tube diameter of 6 mm.



INF150LL

Infusionssystem für die Gabe von Beutelinfusionen in 150 cm Länge mit Luer-Lock Anschlüssen.




UVL

Universal cap for Luer Lock connectors. Available in red.



IKLL-ISO

Latex free injection cap with self-sealing injection membrane for Luer Lock connectors. Available in yellow.



CON-FF1

Circulation connector with 5 cm length with Luer Lock connectors female/female.



CON-MM1

Circulation connector with 5 cm length with Luer Lock connectors male/male.



EXT150LL

Extension line with Luer Lock connectors, 150 cm.



SY100LL

Disposable syringe 100 ml with Luer Lock connector.

Order code	Description	Packaging unit	Connector material
AV-EXT2	Extension lines for double lumen catheters	50 units / carton	PC, PP
AV-EXT2-Y30	Extension lines for dialysis cannulas	80 units / carton	PC, PP
CATH-TAS	DIACATH Catheter dressing	30 units / carton	
CLAMP-SP-02	Stop& Protect compression strap, disposable	50 units / carton	
CON-FF1	Circulation connector, Luer Lock f/f, 5 cm	100 units / carton	PC
CON-MM1	Circulation connector, Luer Lock m/m, 5 cm	100 units / carton	PC
EXT150LL	Extension line with Luer-Lock, 150 cm	25 units / carton	PC
HEDclip	Zugentlastungssystem HEDclip	10 units / carton	
IKLL-ISO	Injection cap LL, sterile, latex-free membrane	100 units / carton	PP
INF150LL	Infusion system with Luer Lock, 150 cm	25 units / carton	PC
KAB75	Sharps container 7,5 l	40 units / carton	
KPA06006-4	Connection and disconnection set for single lumen-catheter (customized sets on request)	20 units / carton	
KPA06004-4	Connection and disconnection set for double lumen-catheter (customized sets on request)	15 units / carton	
KPA06021-6	Connection and disconnection set for shunts	120 units / carton	
K2565	Plastic scissor clamp	50 units / carton	
SNA3	Single-Needle-Adapter with Luer Lock and clamp	50 units / carton	PP, PVC
SY100LL	Disposable Syringe 100 ml	25 units / carton	PP
UVL	Universal cap, Luer Lock, rot	100 units / carton	ABS

